



Gas Law's

Outline

- Introduction
- Charles' Law
- Boyle's Law



Introduction

- **Gas Law's**

There are three gas laws, let start with ideal gas equation.

$$PV = nRT$$

P is pressure of gas

V volume of gas

n is number of moles

T is temperature & **R** is universal gas constant its value is **8.3145J/mol.**

Charles' Law


- Introduction to Charles' Law.

Charles' law states that at a constant pressure, the volume of a fixed mass varies directly with the absolute temperature.

$$V = \text{Const. } T$$

$P = \text{constant}$

- When T increase V increased ($T \uparrow, V \uparrow$)
- When T reduce V reduced ($T \downarrow, V \downarrow$)



Experimental setup

- Equipment Required

- 1) Gas Law Apparatus

- I. Base apparatus


- II. air chamber

- III. Tube for connection

- 1) container of hot water



- 2) container of ice water

- 3) Thermometer



Verification of Charles law.

- 1) Position base apparatus in horizontal direction.
- 2) Connect tube with base apparatus and then to air chamber.
- 3) Put air chamber into hot water.
- 4) Note down position or height of piston (h) and temperature (T).
- 5) Add cold water or ice if available into hot water to lower its temperature.
- 6) Note new temperature and height of piston.
- 7) Repeat **step 5** until you get 6 to 7 values of (T & h) to make graph.



Verification of Charles law.


- $V = \text{Const. } T$
- we need values for Temperature and their corresponding volume.
- We will note down values of compression or expansion to find V by using formula (3.14 r²h) piston.

Diameter of the piston = 32.5 mm


Serial number	Extension <i>mm</i>	Volume <i>mm³</i>	Temperature <i>K</i>
1.	18.3	15204	348
2.	16.1	13390	340
3.	14.5	12028	333
4.	13	10784	327
5.	11.7	9680	322
6.	9.5	7880	313

Verification of Charles law.

- Experimental Results



- Theoretical Results



Verification of Boyle's Law


Introduction to Boyle's Law.

- Boyle's law states that the product of the volume of a gas times its pressure is a constant at a fixed temperature, pressure will be inversely related to the volume

$$V = \text{Const. } 1 / P$$

$T = \text{constant}$

- When V increase P reduced ($P \downarrow$, $V \uparrow$)
- When V reduce P increased ($P \uparrow$, $V \downarrow$)




Apparatus

- Syringe
- Pipe for connection
- Pressure sensor
- Interface
- Computer with data studio




Verification of Boyle's law.

- 1) Setup will be arranged by lab assistant according to figure.
- 2) Open "Data Studio" on computer to which interface is connected.
- 3) Click on create experiment.
- 4) Click on '**start**' button on screen, screen will show graph for **P** and **t**.
- 5) As plunger of syringe pushed pressure increased.
- 6) Note value of volume from syringe and pressure from computer screen.
- 7) Repeat **step 5 to 6** until you get 6 to 7 values of (V & P) to make graph.




Results

- Experimental




- Theoretical



Results

- Graph shown on right is plotted for inverse of Volume.
- There is direct relation between P and $1/V$



Assignment

- Why we have placed base apparatus in horizontal position in charle's law experiment?
- Define 3rd gas law "Gay-Lussac's law"?

Dead line : 22-04-2020 Time 2 O'clock

Email : sperveen.msphy17sns@student.nust.edu.pk

Note

(Do not forget to mention in subject section while composing email)

Group in lab

School & section

Experiment name

